




**valleynorthern**

*Hassle free because we care*

## How to Use and Washing Instructions for Precision Tablet Crusher


### How to Use



*Put pill into Crusher and turn*


### Washing by hand

**1**




*Open Pill Crusher and rinse under tap*

**2**



*Carefully dry with a tissue*

**3**




### Sterilisation

*Steam Time and Temperature*

*18 minutes at 121°C or 3-4 minutes at 134°C*

**1**




**2**




**3**



**4**



**5**



This device complies with ED 1207/2001 and is Year 2000 and Euro compliant